

PowerCure PRO

UPDATED AND STRONGER SMART
HANDHELD LASER DEVICE



Totally Physical
Therapy



Long Duration
Time



Safety & Fast
Efficiency



LCD Sreen




non-invasive

Enjoy the new way of science and technology

650nm+808nm Cold Laser Therapy
Device For Human And Veterinary Use

- ✓ Simplicity of operator
- ✓ Use it twice daily for 15-20 min each
- ✓ Significantly effective






Your physical therapist around you

The curative effect is obvious, with no side effects, small size, light weight, easy to carry, easy to use, and painless treatment

Laser physiotherapy device for pain relief anywhere in the body that you can use at home



Skin Wound Healing



Diminsh Inflammation



Pain Relief




Knee Arthritis



Rheumatoid Arthritis



Sports Injuries & Sprain



Wound & Ulcer
Acupuncture



Soft Tissue Injuries

THERAPY
Application

HOW CAN I GET RELIEF FROM PAIN AT HOME?





- Arthritis/Joint Pain
- Bursitis
- Fibromyalgia
- Temporo-Mandibular Joint Disorder
- Knee Pain
- Finger & Toe Joint Pain
- Muscle Pain /Spasms
- Shoulder /Rotator Cuff
- Tendonitis
- Whiplash
- Back Pain/ Neck Pain
- Carpal Tunnel Syndrome
- Heel Spurs/Plantar Fasciitis
- Migraine Headaches
- Nerve Pain/ Radiculopathy Sciatica
- Sprain/ Strain
- Tennis Elbow
- Trigeminal Neuralgia

A BETTER LIFE FOR YOUR PET

Treat your pet's medical conditions and injuries using PowerCure Laser Therapy Device.









Essential home physiotherapy device for health care

Multiple parts of the whole body apply

Systemic physiotherapy:

Systemic treatment of the joints and acupuncture point, Safe, non-invasive technology without any side effects

Safe and effective:

FDA Cleared Class II Medical Device featuring Laser Therapy, Clinical evidence and safe treatment.

OEM/ODM:

Logo customization acceptable

Important Remarks:

1. Pls treat the injured area at fixed time intervals and directly work on the skin, instead of clothes or bandages.
2. Pls treat the injured area for up to 20-30 minutes, 10 times consecutively, with 5 minute breaks among the treatments.
3. According to the level of pain, continue to treat that area until the problem is solved or an extended relief in pain is experienced, without placing a time limit for each treatment or limiting the treatment period.


Why choose 650nm + 808nm?

650nm Laser penetrates skin to 2~3cm , to relief shallow pain,
808nm Laser can Penetrate to 5~6cm deep area, to treat arthritis
and shoulder pain.



HOW DOES IT WORK?

The product uses 10pcs 650nm wavelength + 5pcs 808nm wavelength laser to illuminate the pain point, promote blood circulation to stimulate cell activity, make the body in a positive state. It produces a range of physiological changes such as improving blood circulation, reducing inflammation, desensitization and pain, wound healing and adjusting body function.



INNOVATION ADVANTAGE

Innovative design is suitable for the whole family

Innovative technology

Innovative spherical lens
New analgesic laser
treatment equipment

- Body sensing comfort
- Unilliating light
- Easy to use




Long battery life

A built-in rechargeable 5000 mah lithium battery, Charge it once and use it for a week

Significant effect

Medical 5pcs 808nm and home 10pcs 650nm laser diode are perfectly combined to effectively treat acute and chronic pain





Total physical therapy

Low-intensity laser therapy technology is a completely physical, safe, non-invasive, no side effects of the technology Use it twice daily for 6-8 min each obvious effect

New magnetic suction structure

Detachable magnetic massage lens, easy to assemble and assemble, convenient to replace other functional treatment head; Direct contact with the skin when in use, easy to operate.






Large LCD screen

The latest pain treatment laser equipment, equipped with ultra big digital screen, More suitable for the elderly

Dual treatment mode on comfort

Continuous mode + pulse mode
Mild release of energy with instantaneous release of stronger energy, penetrate deep muscles, relieve nerves and blood vessels.



SPECIFICATIONS

Laser Medium	GaAlAs semiconductor
Laser Wavelength	808nm \pm 20nm (5PCS) 650nm \pm 20nm (10PCS)
Maximum Output Power per 808nm Laser Diode	250mW \pm 20%
Maximum Output Power per 650nm Laser Diode	5mW \pm 20%
Total Output Power	1300mW \pm 20%
Working Mode	Continuous/Pulse
Time Setting	15 to 60 minutes, 5 minutes as an interval
Battery Capacity	5000mAh lithium battery
Charging Port	DC5V, 2A
Environment Temperature	5 $^{\circ}$ C ~ 40 $^{\circ}$ C
Relative Humidity	20 ~ 80%
Atmospheric Pressure	860hpa ~ 1060 hpa
Adapter Output	DC5V $\overline{=}$ 2A
Adapter Voltage	AC100-240V
Net Weight	270 \pm 5g

GENUINE SECURITY

