


# **SGROW<sup>®</sup> RLT**



- Full Body Red Light Panel
- 660nm 850nm Red Light Therapy



# SGrow 1200 Red Light Panel Specs

Item	1200W Red Light Therapy Panel
Brand Name	SGROW
Power	1200W
Wavelength	660nm:850nm=1:1
Light Color	Red + Infrared
LED QTY	224pcs Dual 5W LED chip
Size	269*874*65mm
Voltage	AC100-277V
Flicker	No Flicker / Flicker-free
Time Control	Yes
LED Control	Yes
Function	Anti aging, pain relief, skin rejuvenation. etc. Health and skin care


## Included:




# ***The Difference of Red Light and Near Infraredlight***

# 660nm Red Light Penetrate to Dermis Layer


# 850nm Infrared Light can deep into Subcutaneous and Muscle Tissue







# 2 Functional Modes VS 4 Functional Mode



“Our Product”


- ❶ Light Mode
- ❷ Timing Mode
- ❸ Environmental Model
- ❹ Recover+ Mode



“Other ”


- ❶ Light Mode
- ❷ Timing Mode

## Multifunctional intelligent panel




- 1** Recover+ Mode:10hz,20hz,40hz.
- 2** Environmental Model: 100%,75%, 50%,25%.
- 3** Timing Mode:20minutes,15minutes, 10minutes,5minutes.
- 4** Light Mode:Red Light/Infrared Light/Red+INF Light

Narrow Beam 30 Degrees  
VS  
Wide Angle Beam 60 Degrees



Narrow Beam 10 Degrees



LED PENAL

Narrow Beam 60 Degrees



LED PENAL

The 60 degree lens has a wider irradiation area than the 10 degree lens, making the irradiation intensity more uniform.

# Hand Hold Design

