









### Say Hello to Your Red-light Therapy Pen!

Each Led Diode is a dual wavelength of 660nm and 850nm. There are three ideal power levels 180mw/c2, 220mw/c2 and pulse mode at 100mw/cm2 for people with sensitive skin.

Note the 850nm is invisible to the human eye.

#### Benefits of 660nm and 850nm

660nm and 850nm is the best wavelength that can penetrate skin and tissue. Clinical trials prove they are the best to relieve pain, help with inflammation, and tissue repair! Take off the focus pen and use it on hands, feet, face, neck, wrists, elbows, knees, ankles, and skin. This device is safe for pets too.

#### Painless Use cases included:

- Speed up Canker Sore Relief
- Face Skin Rejuvenation
- Speed up Cold Sores Relief
- Pain Relief including Teeth, Braces and Mouth
- Inflammation and Open Sores
- Take off the Tip and use the Light on The Entire Body
- Manage Pain and Sores on or around the Genital Area
- Safe for Pets too













